



Random Microseed Matrix-Screening (rMMS), including LCP seeding

Patrick D. Shaw Stewart

Stefan A. Kolek and Patrick D. Shaw Stewart (both Douglas Instruments Ltd, UK), and Bastian Bräuning (Technische Universität München, Germany).

Introduction

Several new techniques that increase the production of crystal structures have emerged in recent years. Random Microseed Matrix-Screening (rMMS), where seed crystals are added automatically to random crystallization screens, is probably the most important of these [1]. During the eight years since the method was published, theoretical understanding of the method has increased [2 - 4], and several important practical variations of the basic method have gained popularity [5, 6]. We will briefly describe some of these variations, including cross-seeding, and introduce a novel method of making LCP seed stocks by scaling up LCP crystallization conditions [7].

We will also describe a method of generating seed gradients across a plate so that the number of crystals in each LCP bolus can be varied, with a practical example.

[1] D'Arcy, A., Villard, F, and Marsh, M. "An automated microseed matrix-screening method for protein crystallization." *Acta Crystallographica Section D: Biological Crystallography* 63.4 (2007): 550-554.

[2] Shaw Stewart, P. D., Kolek, S. A., Briggs, R. A., Chayen, N. E., & Baldock, P. F. (2011). Random microseeding: a theoretical and practical exploration of seed stability and seeding techniques for successful protein crystallization. *Crystal Growth & Design*, 11(8), 3432-3441.

[3] D'Arcy, A., Bergfors, T., Cowan-Jacob, S. W., & Marsh, M. (2014). Microseed matrix screening for optimization in protein crystallization: what have we learned?. *Acta Crystallographica Section F: Structural Biology Communications*, 70(9), 1117-1126.

[4] Shaw Stewart, P. D. & Mueller-Dieckmann, J. (2014). Automation in biological crystallization. *Acta Crystallographica Section F: Structural Biology Communications*, 70(6), 686-696.

[5] Obmolova, G., Malia, T. J., Teplyakov, A., Sweet, R. W., & Gilliland, G. L. (2014). Protein crystallization with microseed matrix screening: application to human germline antibody Fabs. *Structural Biology and Crystallization Communications*, 70(8).

[6] Abuhammad, A., Lowe, E. D., McDonough, M. A., Shaw Stewart, P. D., Kolek, S. A., Sim, E., & Garman, E. F. (2013). Structure of arylamine N-acetyltransferase from *Mycobacterium tuberculosis* determined by cross-seeding with the homologous protein from *M. marinum*: triumph over adversity. *Acta Crystallographica Section D: Biological Crystallography*, 69(8), 1433-1446.

[7] Kolek SA, Braeuning B, Shaw Stewart PD. A novel microseeding method for the crystallization of membrane proteins in lipidic cubic phase. *Acta Crystallographica Section F: Structural Biology Communications*. 2016 Apr 1;72(4):307-12.

I'm going to try to convince you that:

Random microseeding (random Microseed Matrix Screening, or rMMS) should be part of your normal workflow.

You should use the method as soon as possible!

Comment: I'm not talking about "classical" microseeding where you add seed crystals to conditions that are similar to your hit.


Main Reference:

D'Arcy et al. Acta Cryst. (2007). D63:

1. Add seed crystals to a random screen
2. Suspend crushed crystals in the reservoir solution that gave the hits used ("hit solution")
3. Automate!

To get:

- (1) more hits
- (2) better crystals
- (3) the right number of crystals (e.g. for soaking)



Before we start, comments about robotics:

- Contact dispensing is best for microseeding
- It's very helpful if no seed-stock or protein is wasted
- Optimization:
 - o 2-d grid
 - o Combinatorial script
 - o 7-d multivariate designs

Recommended matrix seeding

volumes:

0.3 µl protein

+ 0.2 µl reservoir solution


+ 0.1 µl seed stock

E.g. 3.0 + 2.0 + 1.0 µl for neutron diffraction !

Method of making the seed stock:


See www.douglas.co.uk/mms.htm or sheet

1. Make a seed stock as soon as your crystals stop growing
2. Break crystals with a probe
3. Think about which high-salt solutions are likely to give salt crystals




4. Take 50 μ l of reservoir solution and add everything from a drop
5. Vortex with a Hampton "Seed Bead"
6. Make a dilution series immediately
7. Freeze

Look after your seeds!



1. Dramatic increase in the number of hits
2. The hits are often in completely different conditions (salt/PEG)


Phase diagram of a protein:




Only the blue conditions would be picked up in a typical screen

The red conditions are picked up in rMMS; these are the best conditions to work with

Typical progression of a project:



Cross-seeding, sometimes with mixtures of seed-stocks:




Microseeding can be used with sitting drop ...

... hanging drop ...


... and microbatch-under-oil

How can we use seeding with LCP crystallization?



The ratio of b/a determines the final conc. of precipitant etc. in the LCP. The absolute size of b determines the timing

New LCP hits picked up:



Special parts required for accurate scaling up:

